



● Eppendorf Tubes



● PCR cooler



● PCR plastic products



● Eppendorf Tips



Consumables

It's your sample –
Eppendorf consumables!

eppendorf

	Page
Eppendorf Tubes	
Quality and purity standards	4-5
Technical specifications	6
Safe-Lock Eppendorf Tubes	7
Protein LoBind Tubes	8
DNA/RNA LoBind Tubes	9
Biopur® tubes	10
PCR clean tubes	11
Standard tube 3810X	12
Standard tube 3810	13
Micropestle	13
Eppendorf® Lid _{Bac}	14
epTIPS	
Quality and purity standards	15
epTIPS Standard	16
epTIPS Reloads, epTIPS Box, epTIPS Set	17-18
epTIPS Racks	19
epTIPS Singles	20
epTIPS Filter	21
PerfectPure tips	
PerfectPure tips	22
PerfectPure C-18 tips	23
Cell culture and/or storage of plastic products	
Deepwell plate 96	24
Rack systems	
Tube rack	25
Tube holder 3821	25
PCR plastic products, Combitips and UVette	
Quality standards	26
twin.tec PCR plates	27
0.2 ml and 0.5 ml PCR tubes	28
Work tray and frame, PCR rack	29
Heat-sealing materials	30
PCR film and PCR foil (self-adhesive)	31
Sealing mat 96	32
Cap strips	33
<i>In situ</i> frames	34
PCR cooler	35
Combitips	36
UVette (disposable cuvette for UV/VIS)	37
Appendix	
Centrifugal stability of Eppendorf Tubes	38
Material stability of Eppendorf plastic products	39-43
References	43

Original Eppendorf Tubes

For more than 40 years, the “Eppendorf tip” has proven itself in day-to-day use in the lab – around the world! The increasing sophistication of analytical techniques over the years mean that today’s Eppendorf Tubes must be of higher quality than before.

Top-quality material

All Eppendorf Tubes are made of polypropylene (PP) which meets international purity criteria – for example, the specifications of the BGVV (German regulation for food stuff) for the use of plastics for food handling [1] and the requirements of the FDA. The dyes used for our consumables are free of organic substances and substances containing heavy metals. Quality control ensures that the tubes are checked several times a day for such things as vapor tightness, lid-opening force, force required to pierce the lid, wall strength and centrifugal stability.

Customized purity standards




In correlation with the advanced development of analytical methods, requirements for the purity of consumables has also increased.

For instance, PCR would require a different level of plastic purity than cell culture. Eppendorf Tubes are manufactured in three different quality classes to provide a customized solution for any application: Eppendorf Quality, PCR clean and Biopur®. These products, bearing the seal of Eppendorf quality, meet all the requirements for reliable and easy-to-use consumables. In addition, the products in the PCR clean class are also all free from human DNA, DNase, RNase and PCR inhibitors. Eppendorf Biopur® products are guaranteed sterile, pyrogen-free, RNase-free, DNA-free and ATP-free.

As well as these high Eppendorf quality standards, we also certify compliance with the European Directive on IvD (98/79/EC) for the products listed below. By meeting this standard, Eppendorf offers medical laboratories, in particular, a high degree of reliability for product safety, that guarantees maximum quality, purity and reproducibility for analyses.

IvD
Continuous quality control
Especially suitable for diagnostic laboratories
<ul style="list-style-type: none"> Pipette tips Tubes Combitips plus

Purity standards

	Eppendorf Quality	Eppendorf PCR clean	Eppendorf Biopur®
	Continuous quality control	Testing of production batches (certified)	Batch testing (certified)
	Function, liquid-tightness, precision	Free of human DNA, DNase, RNase and PCR inhibitors	Sterile, free of pyrogens, RNase, DNA and ATP
	<ul style="list-style-type: none"> epTIPS Standard epTIPS Reloads epTIPS Box epTIPS Set Tubes Combitips® plus 	<ul style="list-style-type: none"> epTIPS Reloads epTIPS Filter (PCR clean and sterile) 0.2/0.5 ml PCR tubes PCR plates Capping aids Reaction tubes 	<ul style="list-style-type: none"> epTIPS Racks epTIPS Singles Tubes Combitips® plus
			
Purity standard			
Requirements			
Defined flow properties ¹	●	●	●
Low wettability	●	●	●
High chemical resistance	●	●	●
High thermal stability	●	●	●
Precise mold	●	●	●
ATP-free			● ²
DNA-free (Ph.Eur./USP)		● ¹	● ²
DNase-free		● ¹	●
Individually blister packed			● ³
Free of human DNA		●	●
Free of PCR inhibitors		● ¹	●
Pyrogen-free			● ²
RNase-free		● ¹	● ²
Sterile			● ²

¹ Certified; certificate included with product

² Certified; batch-specific certificate available at www.eppendorf.com


³ Tips also available in boxes/racks

Technical specifications

All Eppendorf standard reaction tubes meet the following quality standards:

- Easy one-handed, multiple opening and closing
- High centrifugal stability across a broad temperature range (see nomogram)*
- Low wettability
- High chemical resistance [2]*
- Autoclavable when open (121 °C, 20 min.)
- Compatible with micropestle (see p. 13)

*See Appendix for further information.



Test conditions:

3810X tubes filled with diluted NaCl solution (concentration 1.2 mg/ml) were rotated in an Eppendorf Centrifuge 5417 R (fixed-angle rotor FA 45-24-11, centrifugable at up to 25,000 x g). To simulate higher rpms and/or greater forces on the tubes, K₂CO₃ and CsCl solutions of various concentrations were used (between 1.4 and 1.9). The 2.0 ml Safe-Lock tubes display a similar response to centrifugation.


Dr. Diana Ludolfs, Bernhardt Nocht Institute Hamburg:

The seal on Safe-Lock Eppendorf tubes is extremely tight thanks to the special lid-locking mechanism.

Safe-Lock Eppendorf Tubes

Product features

- A small hook on the hinged lid clips around the rim of the test tube and prevents the tube from opening accidentally (e.g. during heating)
- Frosted writing surface
- A needle placed through the thin membrane in the middle of the lid allows aerosol-free removal of biohazardous substances.
- Graduation scale



i Safe-Lock Eppendorf Tubes are also available in Biopur®, PCR clean and LoBind quality.


Ordering information

Article	Order no.	Order no.	Order no.
Safe-Lock Eppendorf Tubes	per 500 pcs.	per 1,000 pcs.	per 1,000 pcs.
	0.5 ml	1.5 ml	2.0 ml
colorless	0030 121.023	0030 120.086	0030 120.094
Assortment of yellow, red, blue, green, colorless	0030 121.708	0030 121.694	0030 121.686
Yellow	0030 121.112	0030 120.159	0030 120.205
Red	0030 121.120	0030 120.167	0030 120.213
Blue	0030 121.139	0030 120.175	0030 120.221
Green	0030 121.147	0030 120.183	0030 120.230
Amber (light protection)	0030 121.155	0030 120.191	0030 120.248


Protein LoBind Tubes

Product features

- Minimum protein loss (less than 3% BSA, 1 µg/ml)
- Free of surface coating, e.g. silicone
- Free of DNA, DNase, RNase and PCR inhibitors (PCR clean)
- Special type of polypropylene manufactured under optimized processing conditions
- High degree of transparency



Results: Protein concentration 1 µg/ml



Test conditions:

With ¹²⁵I-labeled BSA, Tris/HCl, pH 7.5, was dissolved in a solution of 0.1% NaCl and incubated at 20 °C for 24 h in the examined tubes. The solutions were then removed and the tubes filled with the above solution without BSA. The radioactivity of the tubes and the solution were then measured. The tubes were rinsed and the radioactivity of the empty tubes measured and recorded.

Ordering information

Article	Order no.
Safe-Lock, Protein LoBind Tubes, PCR clean	
0.5 ml, 250 pieces	0030 108.094
1.5 ml, 250 pieces	0030 108.116
2.0 ml, 250 pieces	0030 108.132

DNA/RNA LoBind Tubes


Product features

- Maximum recovery of DNA/RNA, less than 1% loss at low concentrations under critical conditions (high salt, short DNA fragments (< 150 bp))
- Free of surface coating, e.g. silicone
- Free of DNA, DNase, RNase and PCR inhibitors
- Special polypropylene type, manufactured under optimized processing conditions
- High degree of transparency



i For further information,
www.myliquidtreasures.com

Results: DNA concentration 0.2 ng/µl



Test conditions:

DNA (130 bp) was labeled with ³²P on the 5' end and dissolved in 2.5 M NaCl/TE buffer, pH 7.5. The solution (0.2 ng/µl) was incubated in tubes at 37 °C for 24 h or at 95 °C for 10 min. Following removal of the solutions, the radioactivity of the empty tubes was measured and recorded.

Ordering information

Article	Order no.
Safe-Lock, DNA/RNA LoBind Tubes, PCR clean	
0.5 ml Safe-Lock, 100 pieces	0030 108.035
1.5 ml Safe-Lock, 100 pieces	0030 108.051
2.0 ml Safe-Lock, 100 pieces	0030 108.078


Biopur[®] tubes

Product features

- Sterile
- Free of pyrogen, RNase, DNA and ATP
- Monitored and certified by an external lab
- Individually blistered consumable product for effective protection against follow-up contamination



- 1 Individual batch/lot certificates can be downloaded at www.eppendorf.de



Ms. Joh-Oh Kim, University of Munich Hospital, Grosshadern:

Thanks to the excellent lockable lids, Eppendorf Tubes are reliable and versatile.

Ordering information

Article	Order no.
Biopur[®] Safe-Lock tubes	
0.5 ml, individually sealed, 50 pieces	0030 121.570
1.5 ml, individually sealed, 100 pieces	0030 121.589
2.0 ml, individually sealed, 100 pieces	0030 121.597

PCR clean tubes

Product features

- Free of human DNA and PCR inhibitors
- Free of DNase and RNase
- Certificate available upon request

- 1 For further information, see www.eppendorf.com/pcrclean



Ordering information



Article	Order no.	Order no.	Order no.
	500 pieces	1,000 pieces	1,000 pieces
	0.5 ml	1.5 ml	2.0 ml
Safe-Lock micro test tubes			
PCR clean			
Colorless	0030 123.301	0030 123.328	0030 123.344
Micro test tubes 3810X			
PCR clean			
Colorless		0030 125.215	

Standard tube 3810X

Product features

- Flat, frosted lid and writing surface on body of tube offers a generous labeling area
- Graduation from 0.1–1.5 ml on body of tube

i Eppendorf micro test tubes 3810X are also available in PCR clean quality.



Ordering information


Article	Order no.
Micro test tube 3810X	
1.5 ml, 1,000 pieces	
Colorless	0030 125.150
Blue	0030 125.177
Yellow	0030 125.207
Green	0030 125.185
Red	0030 125.193

Standard tube 3810 Micropestle

Standard micro test tube 3810

Product features

- Ideal for robots, lid is easy to penetrate
- Lower section in lid, small lid area



Ordering information

Article	Order no.
Tube 3810	
1,000 pieces	0030 102.002

Micropestle

Product features

- For resuspending pellets in 1.5 to 2.0 ml test tubes
- Autoclavable (121 °C, 20 min.)




Ordering information

Article	Order no.
Micropestle , for resuspending pellets in 1.5–2.0 ml test tubes,	
Sets of 10	0030 120.973

Eppendorf Lid_{Bac}

Product features

- Tailor-made preparation of culture samples now possible in ml-scale
- No change in procedure – but much easier to handle (cultivation and extraction in the same tube)
- Suitable for freeze-drying and general drying (incubation in vacuum)
- No sterilization of glass equipment
- Ideal, space-saving combinations with Eppendorf Thermomixer comfort as incubator



i For applications and further information, see www.eppendorf.com/lidbac

Ordering information

Article	Order no.
Eppendorf® Lid _{Bac} , contains 100 Safe-Lock 2.0 ml tubes and 100 membrane lids	0030 099.702

Original Eppendorf tips

A pipette is only as good as its tip

Only when using perfectly manufactured pipette tips can the precision of modern pipettes be fully exploited. With this generation of original Eppendorf pipette tips, you can make even better use of this potential. The epTIPS – **Totally Integrated Pipetting System** – have been perfectly matched to Eppendorf pipettes. This results in consistently reduced tip attachment and ejection forces, while ensuring complete sealing. Of course, epTIPS offer benefits when used with pipettes from other manufacturers.

A production philosophy without compromise

In the serial production of epTIPS, strict tolerances set at the boundary of measurable limits is our ideal quality

standard. Since we operate our own production facilities, we can maintain the same high standards in the selection and processing of plastic materials. Ideal wetting properties, high transparency and specially certified purity levels are the visible expression of this production philosophy. The perfect matching of every tip specifically to every Eppendorf pipette guarantees the maximum precision and reliability you can expect from Eppendorf. Eppendorf pipette tips in combination with Eppendorf pipettes satisfy all the requirements of DIN 12 650-2 as well as the imminent standard EN ISO 8655-2.

i For further information, see www.eptips.com

Technical specifications

Eppendorf purity standards

A constantly growing range of applications demands different purity levels with consistently high quality of all products used. Our product range of consumables features three different purity levels: **Eppendorf Quality**, **Eppendorf Biopur** and, specially designed for the high purity requirements of molecular biology labs, **PCR clean**.



For an overview of Eppendorf purity standards, see page 5.

For a complete list of chemical stability properties, see www.eppendorf.com (Support, Applications) under Application No. 5 “Chemical Stability No. 1”.

epTIPS Standard



Product features

- epTIPS Standard are original Eppendorf high-quality pipette tips at an attractive price
- Packed in resealable bags
- Available in all tip sizes from 10 µl to 10 ml
- 200 µl, 300 µl and 1,000 µl tips are color-coded yellow or blue



Ordering information

Article		Order no.
epTIPS Standard		
0.1–10 µl, 34 mm	2 bags of 500 tips = 1,000 tips	0030 000.811
0.1–20 µl, 40 mm	2 bags of 500 tips = 1,000 tips	0030 000.838
0.5–20 µl L, 46 mm	2 bags of 500 tips = 1,000 tips	0030 000.854
2–200 µl, 53 mm	2 bags of 500 tips = 1,000 tips	0030 000.870
20–300 µl, 55 mm	2 bags of 500 tips = 1,000 tips	0030 000.897
50–1,000 µl, 71 mm	2 bags of 500 tips = 1,000 tips	0030 000.919
50–1,250 µl, 76 mm	2 bags of 500 tips = 1,000 tips	0030 000.935
500–2,500 µl, 115 mm	5 bags of 100 tips = 500 tips	0030 000.951
100–5,000 µl, 120 mm	5 bags of 100 tips = 500 tips	0030 000.978
1–10 ml, 165 mm	2 bags of 100 tips = 200 tips	0030 000.765
1–10 ml L, 243 mm	2 bags of 100 tips = 200 tips	0030 000.781

epTIPS Reloads, Box, Set



Product features

- The practical refill system, with trays sorted according to tip size, is packaged either dual-sided or in a 5 stack form
- Contamination-free transfer of trays to the working box
- System optimized for use with multichannel pipettes
- Tips can be used directly from the refill package
- Drastic reduction in waste compared to disposable racks
- All components 100% recyclable



- i** epTIPS Reloads are available in two purity levels: Eppendorf Quality and PCR clean



Ordering information

See next page.

epTIPS Reloads, Box, Set

Ordering information

Article		Order no.
0.1–10 µl, 34 mm	10 trays of 96 tips = 960 tips	0030 073.363
	Reloads, PCR clean, 10 trays of 96 tips = 960 tips	0030 073.746
	Box, 1 reusable box incl. 96 tips	0030 073.002
	Set, 1 reusable box incl. 5 trays of 96 tips	0030 073.207
0.1–20 µl, 40 mm	10 trays of 96 tips = 960 tips	0030 073.380
	Reloads, PCR clean, 10 trays of 96 tips = 960 tips	0030 073.762
	Box, 1 reusable box incl. 96 tips	0030 073.029
	Set, 1 reusable box incl. 5 trays of 96 tips	0030 073.223
0.5–20 µl L, 46 mm	10 trays of 96 tips = 960 tips	0030 073.401
	Reloads, PCR clean, 10 trays of 96 tips = 960 tips	0030 073.789
	Box, 1 reusable box incl. 96 tips	0030 073.045
	Set, 1 reusable box incl. 5 trays of 96 tips	0030 073.240
2–200 µl, 53 mm	10 trays of 96 tips = 960 tips	0030 073.428
	Reloads, PCR clean, 10 trays of 96 tips = 960 tips	0030 073.800
	Box, 1 reusable box incl. 96 tips	0030 073.061
	Set, 1 reusable box incl. 5 trays of 96 tips	0030 073.266
20–300 µl, 55 mm	10 trays of 96 tips = 960 tips	0030 073.444
	Reloads, PCR clean, 10 trays of 96 tips = 960 tips	0030 073.827
	Box, 1 reusable box incl. 96 tips	0030 073.088
	Set, 1 reusable box incl. 5 trays of 96 tips	0030 073.282
50–1,000 µl, 71 mm	10 trays of 96 tips = 960 tips	0030 073.460
	Reloads, PCR clean, 10 trays of 96 tips = 960 tips	0030 073.843
	Box, 1 reusable box incl. 96 tips	0030 073.100
	Set, 1 reusable box incl. 5 trays of 96 tips	0030 073.304
50–1,250 µl, 76 mm	10 trays of 96 tips = 960 tips	0030 073.487
	Reloads, PCR clean, 10 trays of 96 tips = 960 tips	0030 073.860
	Box, 1 reusable box incl. 96 tips	0030 073.126
	Set, 1 reusable box incl. 5 trays of 96 tips	0030 073.320
500–2,500 µl, 115 mm	10 trays of 48 tips = 480 tips	0030 073.509
	Reloads, PCR clean, 10 trays of 48 tips = 480 tips	0030 073.886
	Box, 1 reusable box incl. 48 tips	0030 073.142
	Set, 1 reusable box incl. 5 trays of 48 tips	0030 073.347
100–5,000 µl, 120 mm	Box, 1 reusable box incl. 24 tips	0030 073.169

epTIPS Racks



Product features

- Eppendorf Biopur pipette tips provide maximum biological purity, i.e. guaranteed sterile, pyrogen-free, RNase-free, DNA-free and ATP-free (see page 5), thus satisfying the most stringent demands in medicine, the pharmaceutical and food industries, molecular biology and cell technology
- Continuous control of each batch by an independent laboratory
- Batch-specific certificates available on the internet at www.eppendorf.com
- Packed in racks of 96, 48 or 24 tips



Ordering information


Article		Order no.
epTIPS Racks, Eppendorf Biopur		
0.1–20 µl, 40 mm	40 mm 5 racks of 96 tips = 480 tips	0030 075.005
2–200 µl, 53 mm	53 mm 5 racks of 96 tips = 480 tips	0030 075.021
20–300 µl, 55 mm	55 mm 5 racks of 96 tips = 480 tips	0030 075.048
50–1,000 µl, 71 mm	71 mm 5 racks of 96 tips = 480 tips	0030 075.064
50–1,250 µl, 76 mm	76 mm 5 racks of 96 tips = 480 tips	0030 075.080
500–2,500 µl, 115 mm	115 mm 5 racks of 48 tips = 240 tips	0030 075.102
1–10 ml, 165 mm	5 racks of 24 tips = 120 tips	0030 075.145

epTIPS Singles



Product features

- Eppendorf Biopur pipette tips, individually wrapped in blister packs
- Batch number and expiration date on each blister pack
- Continuous control of each batch by an independent laboratory
- Batch-specific certificates available on the internet at www.eppendorf.com



Ms. Anja Trapp, University of Munich Hospital, Grosshadern:

Eppendorf pipette tips require less force for attachment and ejection forces.

Ordering information

Article	Order no.
epTIPS Singles, Eppendorf Biopur	
0.1–20 µl, 40 mm	100 tips, individually packaged 0030 010.019
2–200 µl, 53 mm	100 tips, individually packaged 0030 010.035
50–1,000 µl, 71 mm	100 tips, individually packaged 0030 010.051

epTIPS Filter



Product features

- Pipette tips with integrated filter made of hydrophobic polyethylene without “self-sealing” additives
- All filter tips are sterile (to USP, DAB, Ph. Eur. standards) and PCR clean (free of human DNA, DNase, RNase and PCR inhibitors)
- Production batch-specific certificates available upon request
- Packed in racks of 96 tips
- If filter becomes wet as a result of incorrect pipetting, it does not swell up and the sample can easily be recovered



i For further information, see www.eppendorf.com/pcrclean



Ordering information

Article	Order no.
epTIPS Filter, PCR clean and sterile	
0.1–10 µl S, 34 mm	10 racks of 96 tips = 960 tips 0030 077.008
0.1–10 µl M, 40 mm	10 racks of 96 tips = 960 tips 0030 077.024
0.5–10 µl L, 46 mm	10 racks of 96 tips = 960 tips 0030 077.040
2–20 µl, 53 mm	10 racks of 96 tips = 960 tips 0030 077.148
2–100 µl, 53 mm	10 racks of 96 tips = 960 tips 0030 077.067
20–300 µl, 55 mm	10 racks of 96 tips = 960 tips 0030 077.083
50–1,000 µl, 76 mm	10 racks of 96 tips = 960 tips 0030 077.105
1–10 ml L, 243 mm	100 tips, individually packaged 0030 077.164

PerfectPure tips

Eppendorf's PerfectPure tips are pipette tips with a flow-through matrix fixed inside the tip outlets. The matrix can function in a wide variety of ways. For instance, a PerfectPure tip can be filled with chromatographic material, and antibodies or other affined reagents and even enzymes can be linked to the matrix. The range of applications is virtually endless: PerfectPure tips are an innovative, high-quality tool for any areas where microliter quantities of valuable samples are purified, concentrated or processed.

As with all of its products, Eppendorf has focused especially on the quality, high sensitivity, specificities and excellent reproducibility of experimental results in the development and manufacturing of PerfectPure tips. As the first product in this new product range, Eppendorf introduces the PerfectPure C-18 tip for preparing peptide samples for mass spectrometry.




PerfectPure C-18 tips

Product features

- Desalts, purifies and concentrates peptide samples for mass spectrometry analyses
- 10 µl pipette tips, functionalized with a flow-through C-18 matrix
- Superior sensitivity supported by optimized products and protocols
- Very good reproducibility thanks to innovative manufacturing technology (patent pending)

Protocol:

2x	Prewetting	50 % acetone
2x	Equilibration	0.1 % TFA
3-10x	Analyte coupling	0.1 % TFA
2x	Wash	0.1 % TFA
1-5x	Elution	Acetone/matrix



Ordering information

Article	Order no.
PerfectPure C-18 tip, 8 tips	0030 008.421
PerfectPure C-18 tip, 96 tips	0030 008.405
PerfectPure C-18 tip, 960 tips	0030 008.413


Deepwell plates 96, 1.2 ml and 2.2 ml

Product features

- Ideal plate for:
 - Cell culture
 - Reaction volumes of up to 1.2 ml or 2.2 ml
 - Storage of samples and transport
- Round well base facilitates recovery of samples
- Compatible with automated systems
- Sealable with self-adhesive seal, heat-sealing or with autoclavable Deepwell mat 96
- Autoclavable (121 °C, 20 min.)



Eppendorf plates articles are made of clean-production polypropylene under clean-room conditions (Class K 100,000) and have been tested to ensure the absence of contaminants. These products meet SBS recommendations to facilitate automated handling.



Ordering information

Article	Order no.
Deepwell plate 96, 1.2 ml, 50 pieces	0030 127.544
Deepwell mat 96, 1.2 ml, 50 pieces	0030 127.552
Deepwell plate 96, 2.2 ml, 50 pieces	0030 127.560
Deepwell mat 96, 2.2 ml, 50 pieces	0030 127.579

Tube rack (multifunctional) Tube holder 3821

Tube rack (multifunctional)

Product features

- For 24 0.5 ml or 1.5–2.0 ml reaction tubes
- Robust coupling mechanism between the individual racks
- Non-slip stand, stackable
- Autoclavable, UV-resistant, labeling surfaces
- Suitable for use in water baths



Ordering information

Article	Order no.
Tube rack, for 0.5 ml tubes	0030 123.107
Tube rack, for 1.5–2.0 ml tubes	0030 123.115

Tube holder 3821

Product features

- For 20 reaction tubes 1.5–2.0 ml
- Numbered bores
- Attractive price



Ordering information

Article	Order no.
Holder 3821, for Eppendorf tubes 3821	3821 000.008

Quality standards

All Eppendorf PCR plastic products satisfy the purity and quality standard of PCR clean. They are certified free of human DNA, DNase, RNase and PCR inhibitors*. Tubes and wells of PCR plates are made of clean-production polypropylene. Polypropylene gives the wells an inert, non-wetting surface that prevents DNA, RNA and enzymes from binding and

improves recovery even when working with the most minute volumes. This provides a much higher product yield than with polycarbonate equipment and reduces the risk of biological contamination.


*Certificates or testing procedures and further information available upon request.



twin.tec PCR plates

Product features

- Excellent rigidity, minimal bending of skirts
- Reduced well-to-well tolerance
- Extremely thin walls for optimum heat transfer
- Diagonalized corner and alphanumeric orientation system



- i** See the twin.tec brochure for complete information.


Ordering information

Article	Ordering No.
twin.tec PCR plate 96, skirted (wells colorless)	
Colorless, 25 pieces	0030 128.648
Yellow, 25 pieces	0030 128.656
Green, 25 pieces	0030 128.664
Blue, 25 pieces	0030 128.672
Red, 25 pieces	0030 128.680
twin.tec PCR plate 96, skirted (wells black)	
Yellow, 25 pieces (not pictured)	0030 128.800
twin.tec PCR plate 96, semi-skirted (wells colorless)	
Colorless, 25 pieces	0030 128.575
Yellow, 25 pieces	0030 128.583
Green, 25 pieces	0030 128.591
Blue, 25 pieces	0030 128.605
Red, 25 pieces	0030 128.613
twin.tec PCR plate 384, (wells colorless)	
Colorless, 25 pieces	0030 128.508
Yellow, 25 pieces	0030 128.516
Green, 25 pieces	0030 128.524
Blue, 25 pieces	0030 128.532
Red, 25 pieces	0030 128.540

0.2 ml and 0.5 ml PCR tubes

Product features


- Special space-saving lid design
 - For maximum capacity in the universal blocks of Eppendorf Mastercyclers
 - For all standard thermocyclers with 0.5 ml block format
- High transparency, even at bottom of tube
- Tight sealing
- Easy to open without contamination
- For use with all standard thermocyclers with 0.2 ml block format



Space-saving lid geometry enables arrangement in the 96-well micro test plate format



Ms. Sigrid Stärk, University of Munich, Oberschleissheim:
Eppendorf PCR tubes are extremely tight-sealing and are easy to label on the frosted surface.



Ordering information

Article	Order no.
0.5 ml PCR tubes, thin-walled with hinged lid , colorless, 500 pieces	0030 124.502
0.2 ml PCR tubes , colorless, 1,000 pieces	0030 124.332
5-tube strip for 0.2 ml PCR tubes , colorless, pack of 125 (= 625 tubes)	0030 124.340
8-tube strip for 0.2 ml PCR tubes , colorless, pack of 120 (= 960 tubes)	0030 124.359

Work trays and frames PCR rack

Work trays and frames

Product features

- 2-part rack for 0.2 ml PCR tubes, consisting of work tray and frame
- Convenient: All tubes can be directly transferred into a 96-well cycler with the removable work tray
- Autoclavable (121 °C, 20 min.)




Ordering information

Article	Order no.
Work trays for 96 0.2 ml PCR tubes , set of 10	0030 124.235
Frame for work tray , set of 5 (together with tray = rack)	0030 124.243

PCR rack

Product features

- Stand platform for handling and storing tubes and plates
- For 0.2 ml, 0.3 ml, 0.5 ml tubes and 96 PCR plates
- Stackable
- Storage down to -90 °C
- Autoclavable (121 °C, 20 min.)



Ordering information

Article	Order no.
PCR rack , 10 pieces	0030 127.455

Heat-sealing materials

Product features

- Hermetic sealing of PCR plates, ideal for low reaction volumes
- Best protection against evaporation during PCR



Name	Heat-sealing film	Peel-it-lite	Pierce-it-lite
Properties	Optically clear polyester/ polypropylene laminate Extremely stable sealing option – cannot be removed or pierced	Laminated aluminum foil Easily removable	Laminated aluminum foil Easily pierced – even with multi-channel pipettes No glue residue on the pipette tips
Seal integrity	80 °C – 140 °C	-200 °C – 120 °C	-80 °C – 120 °C
Sealing time with Eppendorf Heat Sealer	1–2 sec.	2–4 sec.	2–3 sec.
Weldable materials	Polypropylene	Polypropylene Polyethylene	Polypropylene
Special applications	Colorimetric applications Fluorescence applications incl. real-time PCR Storage of hazardous samples	Storage at extremely low temperatures (-200 °C) Foil can be removed at -80 °C Plate is resealable (after removal of the foil through heat-sealing with new foil)	PCR with water bath cyclers Storage and transport of samples

Ordering information

Article	Order no.
Heat-sealing film, 10 x 10 pcs.	0030 127.650
Peel-it-lite, 100 pcs.	0030 127.668
Pierce-it-lite, 100 pcs.	0030 127.676
Foil stripper	0030 127.641

PCR film and PCR foil (self-adhesive)

Product features


- Effective adhesive sealing prevents evaporation loss during reaction and storage
- Removal of film from the plate with no residue PCR film

PCR film

- Sample can be visually monitored through the transparent film

PCR foil

- Easily pierced
- No sticking of the pipette tip (ideal for automated systems)




Ordering information

Article	Order no.
PCR film, (self-adhesive) 100 pcs.	0030 127.480
PCR foil, (self-adhesive) 100 pcs.	0030 127.471

Sealing mat 96

Product features

- Alternative sealing method for 96-well plates
- Multiple use
- Autoclavable (121 °C, 20 min.)



Ms. Tatjana Nonenmacher, University of Regensburg Hospital: *The Eppendorf twin.tec plates are ideally suited for our robots. The plates are available in various colors, making it very easy to differentiate between various reactions.*


Ordering information

Article	Order no.
Sealing mat 96, 5 pcs.	0030 127.463

Cap strips

Product features

- Easy and rapid sealing of Eppendorf twin.tec plates and other standard PCR plates, ideal interface with capping aids
- Easy to remove using extension at end of strip
- Autoclavable (121 °C, 20 min.)



Ordering information

Article	Order no.
Cap strip, (strips of 8), 25 x 12 pcs.	0030 127.498

In situ frames

Product features

- Gas-tight reaction chamber for:
 - *in situ* PCR
 - *in situ* hybridization
 - Microarraying
- Withstands up to 97 °C
- Fits standard microscope slides
- Frame and coverslip can be removed without leaving a residue





Ordering information

Article	Dimensions (mm)	Order no.
<i>In situ</i> frame, 25 µl, 100 pcs.	17 x 17	0030 127.501
<i>In situ</i> frame, 65 µl, 100 pcs.	24 x 24	0030 127.510
<i>In situ</i> frame, 125 µl, 100 pcs.	25 x 36	0030 127.528
<i>In situ</i> frame, 300 µl, 100 pcs.	27 x 68	0030 127.536

PCR cooler

Product features

- Reliable cooling of sensitive samples to 0 °C
- PCR cooler changes colour when temperature exceeds 7 °C
- Samples prevented from freezing
- Easy to use, space-saving and free of contamination because no need for water
- Universal 96-well format for 0.2 ml and 0.5 ml PCR tubes: Strip, plate and individual tubes
- Available in two colors



Ordering information

Article	Order no.
PCR cooler , 0.2 ml starter set (1 x pink, 1 x blue)	3881 000.015
PCR cooler , 0.2 ml pink	3881 000.023
PCR cooler , 0.2 ml blue	3881 000.031

Combitips

Product features

- Ideal system component for all Multipettes and EDOS
- Large selection of volumes with 9 Combitip plus sizes
- 20 different dispensing volumes per Combitip plus size
- Reliable results guaranteed by positive displacement
- Fine tip of Combitip plus guarantees exact dispensing without additional pipette tip
- Available as standard and Eppendorf Biopur® product



Ordering information

Article	Order no. standard	Order no. Biopur
Eppendorf Combitips® plus	100 pcs.	individually wrapped, 100 pcs.
0.1 ml	0030 069.200	0030 069.404
0.2 ml	0030 069.218	0030 069.412
0.5 ml	0030 069.226	0030 069.420
1.0 ml	0030 069.234	0030 069.439
2.5 ml	0030 069.242	0030 069.447
5.0 ml	0030 069.250	0030 069.455
10 ml	0030 069.269	0030 069.463
25 ml	0030 069.293	0030 069.390
50 ml	0030 069.277	0030 069.471
25 ml adapter standard (1 piece), marine blue	0030 069.528	
50 ml adapter standard (1 piece), anthracite	0030 069.161	
25 ml adapter Eppendorf Biopur (set of 7)		0030 069.498
50 ml adapter Eppendorf Biopur (set of 7)		0030 069.480
Variety pack Combitips® plus , standard, 5 pieces of each size plus one 50 ml and one 25 ml adapter	0030 069.285	
Combilong , (2 pcs.) aspirating aid for removing liquids from volumetric flasks and tall bottles for Combitip plus sizes 5, 25, 50 ml	0030 059.506	
Combitip® rack , box with sliding lid for storing up to 100 Combitips plus 0.1 ml to 10 ml	0030 069.897	




Eppendorf Biopur® Combitips plus, individually blister-packed, are guaranteed sterile and free of pyrogen, RNase, DNA and ATP.

UVette* (disposable cuvette for UV/VIS)

Product features

- Suitable for measuring small volumes of 50 µl or more
- Individually blister-packed for contamination-free work
- Free of DNA/RNase and protein
- Two optical path lengths: 2 mm and 10 mm
- UV and VIS-transparent from 220–1,600 nm
- Volume markings for 500 and 1,000 µl
- Fully transparent plastic with outstanding surface properties
- Optimal filling guaranteed by tapered cuvette base
- Low-lying optical window prevents scratches
- Firm-footed design
- Marking possible on frosted grip
- Ideal for use in the BioPhotometer or in standard commercial spectrophotometers with adapter



*US Patent 6249345

Ordering information

Article	Order no.
UVette® , 80 original Eppendorf disposable cuvettes, individually sealed, for direct use in the BioPhotometer without an adapter	0030 106.300
Adapter¹ , for photometers/spectrophotometers with light center height 8.5 mm	4099 001.009
10 mm	4099 002.005
15 mm	4099 003.001
GeneQuant I/II	4099 004.008
Starter set , 80 UVettes + 1 universal adapter for light center height 15 mm (including GeneQuant), convertible to 8.5 mm	4099 100.007
Cuvette stand	4308 078.006

¹Additional adapters for light center heights > 8.5 mm available upon request.

Centrifugal stability of Eppendorf Tubes

Eppendorf tube type	Temperature [°C]	max. RCF [x g]	Time [min.]
Safe-Lock 0.5 ml	40	27000	20
		25000	90
	20	38000	20
		33000	90
Safe-Lock 1.5 ml	40	27000	20
		25000	90
	20	38000	20
		33000	90
Safe-Lock 2.0 ml	40	27000	20
		25000	90
	20	38000	20
		33000	90
3810X	40	27000	20
		25000	90
	20	38000	20
		33000	90
3810	40	27000	20
		25000	90
	20	38000	20
		33000	90

Test conditions:

3810X tubes filled with diluted NaCl solution (concentration 1.2 mg/ml) were rotated in an Eppendorf Centrifuge 5417 R (fixed-angle rotor FA 45-24-11, centrifugable at up to 25,000 x g). To simulate higher rpms and/or greater forces on the tubes, K₂CO₃ and CsCl solutions of various concentrations were used (between 1.4 and 1.9). The 2.0 ml Safe-Lock tubes display a similar response to centrifugation.

Please note that high-density liquids (e.g. chloroform) reduce the max. time or max. RCF at which the tubes can be centrifuged. Organic solvents and acids may reduce the stability of the plastic. For more details, please refer to the accompanying material stability overview.

The degree of mechanical load to which Biopur tubes can be subjected during centrifugation is reduced. If in doubt, please contact Eppendorf.

Tip for ultracentrifugation: It is possible to centrifuge up to 70,000 x g when form-fitting suspension racks are used and the tubes (except PCR tubes) are filled to their nominal volume.

Technical information subject to change.

Material stability of Eppendorf plastic products

1 = stable; life of several months

2 = conditionally stable; life of a few weeks

3 = unstable; life of a few hours and/or rapid disintegration

0 = not tested

Chemical	Concentration %	PP		
		+20 °C	+40 °C	+60 °C
A				
Acetaldehyde, aqueous solution	40	1	1	1
Acetic acid, aqueous solution	25–60	1	1	1
Acetone (dimethyl ketone)	100	1	1	2
Acrylonitrile		1	0	0
Adipic acid, aqueous solution	saturated	1	1	1
Allyl alcohol (2-propene-1-ol)	96	1	1	1
Aluminum chloride, aqueous solution	saturated	1	1	1
Aminic acid	100	1	1	2
Ammoniac, aqueous solution	saturated	1	1	1
Ammonium chloride, aqueous solution	saturated	1	1	1
Ammonium hydroxide, aqueous solution	30	0	0	2
Amyl acetate	100	2	0	3
Amyl alcohol (1-Pentanol)	100	1	1	1
Amyl chloride (1-Chloropentane)	100	3	0	0
Aniline	100	1	1	1
Aniline, aqueous solution	saturated	2	2	2
Aqua regia (HNO ₃ , conc. HCl)		2	0	3
B				
Benzaldehyde, aqueous solution	saturated	1	0	0
Benzene	100	2	0	3
Benzol	100	2	0	3
Benzoyl chloride		2	0	0
Benzyl alcohol	100	1	0	2
Boric acid, aqueous solution	saturated	1	1	1

All information provided represents a recommendation without guarantee.

Chemical	Concentration %	PP		
		+20 °C	+40 °C	+60 °C
B				
Bromic acid		1	0	0
Bromine, liquid	100	3	3	3
Bromobenzene		3	3	3
1,3-butadiene	100	2	0	3
Butanoic acid	100	1	0	0
Butanoic acid, aqueous solution	20	1	0	0
Butylacetate	100	2	0	3
C				
Calcium chloride, aqueous solution	saturated	1	1	1
Calcium hydroxide, aqueous solution	any	1	1	1
Calcium hypochloride, 12.5% active chlorine		1	1	1
Carbon disulphide	100	1	0	3
Carbon tetrachloride (tetra-chloromethane)	100	3	3	3
Chlorobenzene	100	1	0	0
Chloroform (trichloromethane)	100	2	0	3
Chromic acid, aqueous solution	50	2	2	2
Chromic/sulfuric acid		3	3	3
Citric acid	any	1	1	1
Copper sulfate, aqueous solution	saturated	1	1	1
Cresol, aqueous solution	up to 90	1	0	0
Cyclohexane	100	1	1	1
Cyclohexanol	100	1	0	2
Cyclohexanone	100	1	0	2
D				
Decahydronaphthalene	100	2	2	2
Dibutyl phthalate	100	1	0	2
Dichlorobenzene		2	2	0
Diesel fuel		1	0	2
Diethylene glycol		1	1	1
Diethylether		2	0	0
Diisopropylether	100	2	0	3
Dimethylformamide	100	1	1	1
Di-N-butyl ether		2	0	3
Dioxane	100	2	2	2

All information provided represents a recommendation without guarantee.

Chemical	Concentration %	PP		
		+20 °C	+40 °C	+60 °C
E				
Ethanol, aqueous solution	any	1	1	1
Ethanol	100	1	1	1
Ethyl acetate	100	1	2	2
Ethyl benzene	100	2	0	0
Ethylene chloride (1,2-dichloroethane)	100	2	0	3
Ethylene oxide (1,2-epoxy ethane)	100	1-2	0	0
F				
Formaldehyde (methanal), aqueous solution	40	1	1	1
G				
Glycerine	100	1	1	1
Glycol	100	1	1	1
Glycol, aqueous solution	any	1	1	1
H				
Heating oil	100	1	0	2
Heptane	100	2	2	2
Hexane	100	1	0	2
Hydrochloric acid, aqueous solution	over 30	1	1	1
Hydrogen peroxide, aqueous solution	up to 30	1	0	2
Hydrogen peroxide, aqueous solution	90	1	0	0
Hydrofluoric acid	up to 40	1	1	1
Hydroquinone (1,4-dihydroxybenzene)	any	1	1	1
I				
Iodine-iodine-potassium solution	50	1	1	1
Isopropanol, aqueous solution	any	1	1	1
Isopropanol	100	1	1	1
L				
Lactic acid (2-hydroxypropane acid), aqueous solution	up to 90	1	1	1
Linseed oil	100	1	1	1

All information provided represents a recommendation without guarantee.

Chemical	Concentration %	PP		
		+20 °C	+40 °C	+60 °C
M				
Mercury	100	1	1	1
Methanol	100	1	1	1
Methyl acetate	100	1	1	1
Methylene chloride (dichloromethane)	100	2	3	3
Mineral oil	100	1	0	2
Monochloroacetic acid, aqueous solution	85	1	1	1
Monochloroacetic acid	100	1	1	1
N				
Nitric acid, aqueous solution	up to 30	1	0	2
Nitric acid, aqueous solution	65	3	3	3
Nitrobenzene	100	1	1	1
O				
Oleic acid (cis 9-octadecenoic acid)	100	1	0	2
Oxalic acid, aqueous solution	any	1	1	2
P				
Paraffin oil	100	1	0	2
Pechloric acid, aqueous solution	20	1	1	1
Petroleum ether	100	1	0	2
Petroleum		2	0	0
Petroleum	100	1	0	2
Phenol, aqueous solution	90	1	1	1
Phosphoric acid, aqueous solution	80	1	1	1
Phosphoric acid, aqueous solution	95	1	0	2
Potassium chloride, aqueous solution	saturated	1	1	1
Potassium permanganate, aqueous solution	saturated	1	1	1
Pyridine	100	2	2	2
S				
Silicone oil	100	1	1	1
Silver nitrate, aqueous solution	any	1	1	1
Sodium acetate, aqueous solution		1	1	1
Sodium hypochloride, aqueous solution	diluted	1	1	1-2
Sulfuric acid, aqueous solution	50	1	1	1
Sulfuric acid, aqueous solution	96	2	0	3

All information provided represents a recommendation without guarantee.

Chemical	Concentration %	PP		
		+20 °C	+40 °C	+60 °C
T				
Tartaric acid (butanedioic acid, 2,3 – dihydroxy), aqueous solution	saturated	1	1	1
1,1,2,2-tetrachloro-ethane	100	2	0	3
Tetra-hydrofuran	100	2	0	3
U / V				
Urea, aqueous solution	saturated	1	1	1
Vaseline		1	0	1-2
X / Z				
Xylene	100	3	3	3
Zinc chloride, aqueous solution	diluted			
Zinc chloride, aqueous solution	saturated	1	1	1
Zinc sulfate, aqueous solution	diluted	1	1	1
Zinc sulfate, aqueous solution	saturated	1	1	1

All information provided represents a recommendation without guarantee.

Source: [2]

The Eppendorf twin.tec PCR plates are a combination of two materials, polycarbonate (in the frame) and polypropylene (in the wells), which enables the plates to be used in molecular biology applications based on their very good chemical stability against typical PCR reagents and chemicals, e.g. DMSO (10%, for RT), isopropanol and glycerin.

For a complete list on chemical stability, see www.eppendorf.com (Support, Applications) in Application No. 5 “Chemical Stability No. 1”.

References

- [1] GdV Form Umweltrecht; Bedarfsgegenständeverordnung (i.e. “German Consumer Goods Ordinance”); Version of 20 June 2002.
 [2] Eppendorf Application No. 56; 03/03: “The best material for original Eppis!”



Research
Thermomixer
MiniSpin®
Multipipette®
BioPhotometer
Research pro
epMotion 5070
Mastercycler®



eppendorf

In touch with life

Your local distributor: www.eppendorf.com/worldwide

Application Hotline: +49 180 3666789

Eppendorf AG · 22331 Hamburg · Internet: www.eppendorf.de · Phone: +49 (0)40 538 01 0

eppendorf® is a registered trademark. · Printed in Germany
Order no. A003 113 021/GB1/10T/0204/NEUH · Chlorine-free bleached paper