

810-SG Series **COMPACT GLOVE BOX** with transfer chamber AND flat side access door

Completely sealed glove box for basic isolation and/or containment.



The 810-SG Series by PLAS■LABS are constructed of an optically clear material, developed for isolation and / or containment. All units feature a high visibility of the interior, cleanliness, safety, complete containment, portability and a draft-free atmosphere. The installed purging valves are ideal for lowering oxygen and humidity levels.

Other features include:

- Transfer chamber with purging valves and vacuum gauge.
- Flat side access door for easily introducing larger equipment.
- Four (4) key cock valves for purging: Two on the main body, two on the transfer chamber.
- Hospital grade multiple electrical outlet strip.
- Pressure relief valve with small HEPA filter.
- Standard 24 month warranty, (not on gloves or consumables.)
- White ambidextrous gloves.

	810-SG10	810-SG20	810-SG30	810-SG40
Inside Dimensions	28"w x 23"d x 29"h 710mm x 580mm x 740mm	35"w x 29"d x 30"h 890mm x 740mm x 760mm	48"w x 29"d x 32"h 1220mm x 740mm x 810mm	60"w x 29"d x 32"h 1524mm x 740mm x 810mm
Outside Dimensions (Includes transfer chamber 12" long)	43"w x 24"d x 31"h 1100mm x 610mm x 790mm	49"w x 30"d x 31"h 1250mm x 760mm x 790mm	63"w x 31"d x 35"h 1600mm x 790mm x 860mm	76"w x 32"d x 35"h 1930mm x 812mm x 890mm
Access Door Opening	15.5" x 22"h 400mm x 560mm	21.5" x 22"h 546mm x 560mm	21.5" x 24"h 546mm x 610mm	21.5" x 24"h 546mm x 610mm