

Petri Dishes

I. APPLICATIONS

1. Microbiology
2. Microbiological culture and automatic filling

II. FEATURES

1. Materials: medical grade polystyrene
2. EO(Ethylene Oxide) sterilized
3. Higher clarity than glass

III. PETRI DISH DESIGN

1. Available in 90 x 16mm
2. Vented lid design for gas exchange
3. Slideable
4. Accurate dimensions for automatic filling devices

IV. PACKAGING

1. Bulk, 500 pieces/case
2. Sterile

V. DESCRIPTION

Item	Description	Packaging	Sterile
66-1701	90 x 16mm petri dish	Bulk, 10pieces/bag, 50bags/case	Yes

VI. STERILITY

1. Products are EO (Ethylene Oxide) sterilized according to ISO13485 standard.
2. For each lot, sterilization records are recorded and signed by qualified personnel for product release.

VII. CERTIFICATION



Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

VIII. RESTRICTIONS


Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

VII. PRODUCT PARAMETER



Dimensions=mm




SECTION B-B




Lid



SECTION A-A



Dish



Set