

## Racked Filter Tips

### I. APPLICATIONS

1. For use with pipettors for the measure and transfer of liquids/chemicals

### II. FEATURES

1. Materials: Tips are made of PP, filters are made of UHMWPE.
2. Chemically resistant to alcohols and mild organic solvents.
3. 100% inert.
4. Sterile. Certified RNase & DNase Free.
5. Low retention treated tips undergo a coating-free treatment and will not contaminate samples.

### III. FILTER TIPS DESIGN

1. 10 $\mu$ l, 20 $\mu$ l, 100 $\mu$ l, 200 $\mu$ l, and 1000 $\mu$ l tip sizes.
2. Sterile racked packaging.
3. Filters eliminate errors caused by saturation vapor pressure.

### IV. DESCRIPTION

Cat. No.	Description
23-0010A	10 $\mu$ l, rack pack. 96 Pieces/Rack, 100 Racks/Case
23-0010H	10 $\mu$ l, rack pack. Extra-long, 96 Pieces/Rack, 50 Racks/Case
23-0011H	10 $\mu$ l, rack pack, low-retention, 96 Pieces/Rack, 50 Racks/Case
23-0020A	20 $\mu$ l, rack pack, 96 Pieces/Rack, 100 Racks/Case
23-0100	100 $\mu$ l, rack pack. 96 Pieces/Rack, 100 Racks/Case
23-0100A	100 $\mu$ l, rack pack. 96 Pieces/Rack, 100 Racks/Case
23-0200A	200 $\mu$ l, rack pack, 96 Pieces/Rack, 100 Racks/Case
23-0200A1	200 $\mu$ l, rack pack. 96 Pieces/Rack, 50 Racks/Case
23-0200A2	200 $\mu$ l, extra-long, rack pack. 96 Pieces/Rack, 50 Racks/Case
23-0200H	200 $\mu$ l, extra-long, rack pack. 96 Pieces/Rack, 50 Racks/Case
23-0201H	200 $\mu$ l, extra-long, rack pack, low-retention, 96 Pieces/Rack, 50 Racks/Case
23-1000	1000 $\mu$ l, rack pack, 96 Pieces/Rack, 60 Racks/Case
23-1000A	1000 $\mu$ l, extra-long, rack pack, 96 Pieces/Rack, 50 Racks/Case
23-1001A	1000 $\mu$ l, extra-long, rack pack, low retention, 96 Pieces/Rack, 50 Racks/Case
23-1000H	1000 $\mu$ l, rack pack, 96 Pieces/Rack, 50 Racks/Case
23-1250H	1000 $\mu$ l, extra-long, rack pack, 96 Pieces/Rack, 50 Racks/Case
23-1251H	1000 $\mu$ l, extra-long, rack pack, low retention, 96 Pieces/Rack, 50 Racks/Case

### V. STERILITY Yes

### VI. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

## VII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

## VIII. Photos



**23-0010A**



**23-0010H**



**23-0011H**



**23-0020A**



**23-0100**




**23-0100A**



**23-0200A**



**23-0200A1**



**23-0200A2**



**23-0200H**



**23-0201H**



**23-1000**



**23-1000A**



**23-1001A**




**23-1000H**




**IX. PRODUCT PARAMETER**

**Dimensions=mm**




23-0010A




23-0010H  
23-0011H




23-0020A




23-0100




23-0100A



23-0200A




23-0200A1




23-0200A2




23-0200H  
23-0201H




23-1000



23-1000A  
23-1001A



23-1000H



23-1250H  
23-1251H