

## Cell Culture Flasks

### I. APPLICATIONS

- 1.Suitable for cell and tissue culture

### II. FEATURES

- 1.Body: High quality medical grade polystyrene; Cap: High-density polyethylene; Filter: Acrylic copolymer, 0.22µm
- 2.Sterile
- 3.Non-pyrogenic
- 4.Non-cytotoxic
- 5.DNase, RNase, and human DNA free
- 6.Stackable
- 7.Surface Treatment

### III. DESIGN

- 1.Cell culture surface area of 25 cm<sup>2</sup>, 75 cm<sup>2</sup>, or 175 cm<sup>2</sup>
- 2.Plug cap and filter cap options available for all flask sizes
- 3.Short, wide, angled neck design for easy access

### IV. PACKAGING

1. Bulk,5 Pieces/bag, 200 Pieces/case.(07-8025,07-9025)
2. Bulk,5 Pieces/bag, 100 Pieces/case.(07-8075,07-9075)
3. Bulk,5 Pieces/bag, 40 Pieces/case.(07-8175,07-9175)

### V. DESCRIPTION

Item	Cap type	Culture Area (cm <sup>2</sup> )	Working Volume (ml)	Total Volume (ml)	Packaging
07-8025	Filter	25	7	60	5 Pieces/bag, 200 Pieces/case
07-8075	Filter	75	25	250	5 Pieces/bag, 100 Pieces/case
07-8175	Filter	175	50	650	5 Pieces/bag, 40 Pieces/case
07-9025	plugged	25	7	60	5 Pieces/bag, 200 Pieces/case
07-9075	plugged	75	25	250	5 Pieces/bag, 100 Pieces/case
07-9175	plugged	175	50	650	5 Pieces/bag, 40 Pieces/case

### VI. STERILITY

1. Sterile

### VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard

**VIII. RESTRICTIONS**

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

**IX. PRODUCT PARAMETER**

Dimensions=mm

Cat #	Description	A	B	C	D	E
<b>07-8025</b>	Cell Culture Flask,Filter Cap	50.10	72.00	92.38	14.50	24.60
<b>07-8075</b>	Cell Culture Flask,Filter Cap	84.20	124.10	151.04	19.00	32.00
<b>07-8175</b>	Cell Culture Flask,Filter Cap	116.60	199.08	226.69	20.50	34.72
<b>07-9025</b>	Cell Culture Flask,Plug Cap	50.10	72.00	92.38	14.50	24.60
<b>07-9075</b>	Cell Culture Flask,Plug Cap	84.20	124.10	151.04	19.00	32.00
<b>07-9175</b>	Cell Culture Flask,Plug Cap	116.60	199.08	226.69	20.50	34.72

