

## Cell Scrapers & Lifters

### I. APPLICATIONS

1. Cell culture
2. Cutting and scraping of surfaces

### II. FEATURES

1. Materials: polyethylene
2. Gamma radiation sterilized
3. Non-pyrogenic

### III. SCRAPER/LIFTER DESIGN

1. Scrapers available in 18cm and 25cm, lifter is 18cm
2. 1.8cm blade length for scraper, 2.0cm blade length for lifter
3. Pivoting blade on scraper for precise cutting and scraping
4. Individually wrapped

### IV. PACKAGING

1. Bulk, 100 pieces/case or 200 pieces/case

### V. DESCRIPTION

Item	Description	Packaging	Sterile
70-1180	18cm cell scraper	Bulk, 100 pieces/case	Yes
70-1250	25cm cell scraper	Bulk, 200 pieces/case	Yes
70-2180	18cm cell lifter	Bulk, 100 pieces/case	Yes

### VI. STERILITY

1. Products are Gamma Radiation Sterilized according to ISO13485 standard.
2. For each lot, sterilization records are recorded and signed off by qualified personnel for product release.

### VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

### VIII. RESTRICTIONS

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

### IX. PRODUCT PARAMETER

**Dimensions=mm**

<b>Cat #</b>	<b>Description</b>	<b>A</b>	<b>B</b>	<b>C</b>	<b>D</b>	<b>E</b>
<b>70-1250</b>	Cell Scraper	250.0	18.0	9	89.0	O.D. 7.7
<b>70-1180</b>	Cell Scraper	180.0	18.0	9	85.0	O.D. 7.4
<b>70-2180</b>	Cell Lifter	180.0	20.0	7.4	-	-

