

Haier Biomedical

Intelligent Protection of Life Science

-25°C Biomedical Freezer

DW-25L262



Ergonomic Design



Lockable design



Intelligent alarm



Reliable Performance

High density foam insulation for rigidity and stable storage temperature
LED digital display for clear observation



Adjustable Temperature

Adjustable freezer compartment range:
-10 °C ~ -25 °C



Two Alarm Functions

High/low temperature, sensor error

Specifications

DW-25L262

Technical Data

Type	-25°C Deep Freezer
Cabinet Type	Upright
Climate Class	4
Cooling Type	Direct cooling
Defrost Mode	Manual
Refrigerant	CFC-Free
Noise Level (dB(A))	34

Performance

Cooling performance(°C)	-25
Temperature Range(°C)	-10~-25

Control

Controller	Microprocessor
Display	LED

Electrical Data

Power Supply(V/Hz)	220~240/50
Power(W)	90
Electrical Current(A)	0.7

Dimensions

Capacity(L/Cu.Ft)	262/9.3
Net/Gross Weight(approx)	88/93(kg) 194.0/205.0(lbs)
Interior Dimensions(W*D*H)	480*465*1430(mm)18.9*18.3*56.3(in)
Exterior Dimensions(W*D*H)	700*705*1665(mm) 27.6*27.8*65.6(in)
Packing Dimensions(W*D*H)	775*760*1880(mm) 30.5*29.9*74.0(in)
Container load (20'/40'/40'H)	24/48/48



Functions

High/Low Temperature	Y
Sensor Error	Y

Others

Certificate	CE
-------------	----

Haier Medical and Laboratory Products Co., Ltd

Brand Building, Haier Industrial Park, No. 1
Haier Road, Qingdao, 266101, P.R. China
Tel: +86-0532-88936011 / 5955
Website: www.haiermedical.com



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.

