

# Haier Biomedical

Intelligent Protection of Life Science

## Combined Refrigerator and Freezer

HYCD-282



### Refrigerator & Freezer Combination

Independent display for refrigerator section and freezer section



### Alarm Functions

High/low temperature, power failure, sensor error, low battery, high/low ambient temperature, and remote alarm



### Standard Stainless Steel Interior

Both refrigerator and freezer chamber are stainless steel interior, durable, easy to clean and anti-corrosion

### Ergonomic Design



Lockable design



Intelligent alarm



Foot/Caster design



Test hole design

# Specifications

## HYCD-282

### Technical Data

Type	Combined Refrigerator and Freezer
Cabinet Type	Upright
Climate Class	4
Cooling Type	Refrigerator: Forced air cooling Freezer: Direct cooling
Defrost Mode	Refrigerator: Auto Freezer: Manual
Refrigerant	CFC-Free
Noise Level (dB(A))	43

### Performance

Temperature Range(°C)	Refrigerator: 2~8 Freezer: -20~-40
-----------------------	---------------------------------------

### Control

Controller	Microprocessor
Display	LED

### Electrical Data

Power Supply(V/Hz)	220-240/50
Power(W)	400
Electrical Current(A)	2.8

### Dimensions

Capacity(L/Cu.Ft)	Refrigerator:185/6.5 Freezer:97/3.4	
Net/Gross Weight(approx)	145/160(kg)	319.7/352.7(lbs)
Interior Dimensions(W*D*H)	Refrigerator:605*510*720(mm) Freezer:515*465*440(mm)	23.8*20.0*28.3(in) 20.3*18.3*17.3(in)
Exterior Dimensions(W*D*H)	736*660*1810(mm)	29.0*26.0*71.3(in)
Packing Dimensions(W*D*H)	775*720*1930(mm)	30.5*28.3*76(in)
Container load (20'/40'/40'H)	21/45/45	



### Functions

High/Low Temperature	Y
Remote Alarm	Y
Power Failure	Y
Sensor Error	Y
Low Battery	Y
High Ambient Temp	Y

### Accessories

Caster	Y
Foot	Y
Test Hole	Y
Shelves/Drawers	Refrigerator: 3/1 Freezer: 0/2
USB Interface	Y

### Others

Certificate	CE
-------------	----

## Haier Biomedical UK

Address: 2 Birch House, Brambleside, Bellbrook  
Industrial Estate, Uckfield, East Sussex,  
TN22 1QQ, United Kingdom  
TEL: +44(0) 1825 744 696 (Head Office)  
[www.haierbiomedical.co.uk](http://www.haierbiomedical.co.uk)



Note: If slight difference occurs between pictures and actual products, please refer to actual products. Our company reserves the right of final interpretation to this brochure and please contact with us for further information you might need.



Scan 2D barcode with your mobile device to learn more information about Haier Biomedical