

Haier

R1B1 Infrared Thermometer Introduction


海尔集团
中国家电第一品牌

R1B1-- Product Features

Haier



- ✓ Measure in 1s
- ✓ Auto power-off




- ✓ Memory recall
- ✓ Fever warning

- ✓ Non-touch measurement
- ✓ Forehead mode



- ✓ Practical patent design
- ✓ High definition display



R1B1--Product Parameters

Haier



ITEM	R1B1
Product Size	37*37*138mm
Product Weight	47g
Memory	32 groups
Operating temp	5°C~40°C
Measuring mode	Forehead mode
Measurement range	32.0°C~42.2°C
Accuracy	±0.2°C(35.0°C~42.0°C) ±0.3°C(32.0°C~34.9°C /42.1°C~42.2°C)
Voltage source	2*1.5V Batteries size AAA
Automatic power off	60s

R1B1-- Product Packaging

Haier






ITEM	R1B1
BOX SIZE	65*45*190mm
CTN SIZE	470*410*405mm
CTN	120pcs
CMB	0.078m ²
G.W	10.4KG



R1B1 Qualification—CE/FDA

Haier

EC Certificate

Production Quality Assurance System
Directive 93/42/EEC on Medical Devices (MDD), Annex V
(Devices in Class IIa, IIb or II)

No. G2 078453 0014 Rev. 01

Manufacturer: Alicn Medical Shenzhen, Inc
4/F, B Building, Shenfubao Modern Optical Factory, Kengzi Street, Pingshan District, 518122 Shenzhen City, PEOPLE'S REPUBLIC OF CHINA

Facility(ies): Alicn Medical Shenzhen, Inc
4/F, B Building, Shenfubao Modern Optical Factory, Kengzi Street, Pingshan District, 518122 Shenzhen City, PEOPLE'S REPUBLIC OF CHINA


Product Category(ies): Clinical Infrared Thermometer, Electronic Thermometer and Blood Pressure Monitor

The Certification Body of TÜV SÜD Product Service GmbH declares that the aforementioned manufacturer has implemented a quality assurance system for manufacture and final inspection of the respective devices / device categories in accordance with MDD Annex V. This quality assurance system conforms to the requirements of this Directive and is subject to periodical surveillance. For marketing of class IIb and III devices an additional Annex III certificate is mandatory. See also notes overleaf.


Report No.: G21805201

Valid from: 2019-07-15
Valid until: 2022-05-07

Date: 2019-07-15


 Stefan Probst
 Head of Certification/Notified Body

Page 1 of 1
TÜV SÜD Product Service GmbH is Notified Body with identification no. 0123
TÜV SÜD Product Service GmbH - Certification Body - Hollenstraße 65 - 80339 Munich - Germany



U.S. FOOD & DRUG ADMINISTRATION

[Follow FDA](#) | [En Español](#)

SEARCH

Home
Food
Drugs
Medical Devices
Radiation-Emitting Products
Vaccines, Blood & Biologics
Animal & Veterinary
Cosmetics
Tobacco Products

Establishment Registration & Device Listing

[FDA Home](#)
[Medical Devices](#)
[Databases](#)

1 result found for Establishment Registration or Business Trade Name : Alicn New Search

Establishment Name	Registration Number	Current Registration Yr
ALICN MEDICAL SHENZHEN, INC CHINA	3009889471	2019
<ul style="list-style-type: none"> • Laryngoscope, Rigid Manufacturer • Thermometer, Electronic, Clinical - Non-Contact Infrared Thermometer AET-R161, AET-R171, AET-R181, AET-R1D2 And AET-R1D1; Non-Contact Infrared Thermometer AET-R1B6 Contract Manufacturer, Manufacturer • System, Measurement, Blood-Pressure, Non-Invasive - Arm Blood Pressure Monitor Models Of AES-U111, AES-U141, AES-U511, AES-U312, AES-U512, AES-U531, Arm Blood Pressure Monitor, Models Of AES-U112, AES-U122, AES-U151, AES-U522, AES-U311 And AES-U321; Arm Blood Pressure Monitor, Models Of AES-U121, AES-U331, AES-U341, AES-U351, AES-U361, AES-U512; Arm Blood Pressure Monitor, Models Of AES-U171, AES-U181, AES-U131; Arm Blood Pressure Monitor, Models Of AES-U631, AES-U632, AES-U532, AES-U132, AES-U612, AES-U611 Contract Manufacturer, Manufacturer 		

Can't find what you're looking for? [Try a new search](#)

Page Last Updated: 12/02/2019
 Note: If you need help accessing information in different file formats, see Instructions for Downloading Viewers and Players.
 Language Assistance Available: Español | 繁體中文 | Tiếng Việt | 한국어 | Tagalog | Русский | العربية | Kreyòl Ayisyen | Français | Polski | Português | Italiano | Deutsch | 日本語 | فارسی | English

Haier

Thank You

海尔集团
中国定制 专用指定文件