

Microscope Slides

I. APPLICATIONS

1. Suitable for application in automatic machinery.

II. FEATURES

1. The slides are with clean glass surface
2. Ground edges prevent personal injury and chipping of glass caused by the cut edges
3. The glass slide end faces are polished
4. The slides without frosted area can be effectively used
5. The slides with frosted area are convenient for classification and storage

III. SLIDES DESIGN

1. Size: 25.4 x 76.2mm
2. Two kinds: Plain Slides and Single Frosted End
3. Materials: Clear glass

IV. PACKAGING

1. Bulk, 72 Pieces/pack, 50 Packs/case

V. DESCRIPTION

Item	Description	Packaging	Sterile
7101	Plain Slides	72 Pieces/pack, 50 Packs/case	No
7105	Single Frosted End	72 Pieces/pack, 50 Packs/case	No

VI. STERILITY

1. N/A

VII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

VIII. RESTRICTIONS:

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

IX. PRODUCT PARAMETER

Dimensions=mm

Cat #	Description	A	B	C	D	E
7101	Plain Slides	25.4	76	2	25.4	-
7105	Single Frosted End	25.4	76	2	25.4	20

