Document Number: **BD-T0060** 

Document Date: May 2012

Page 1/3

# **PCR Strip Tubes**

#### I. APPLICATIONS

- 1. Molecular biology
- 2. Cell centrifugation
- 3. Sample storage
- 4. Centrifugation of precipitates

## **II. FEATURES**

- 1. Materials: Virgin polypropylene
- 2. Chemically resistant to alcohols and mild organic solvents
- 3. Transparent wall permits easy viewing of tube contents
- 4. Temperature range: stable from -20°C to 100°C
- 5. Tested to withstand centrifugation of 6,000 RCF in a fully supported rotor
- 6. Certified RNase & DNase Free
- 7. Extra transparent caps allow for RT-PCR applications

#### **III. TUBE DESIGN**

- 1. Available in 0.2ml version
- 2. Conical bottom design fits most standard floor model and table-top centrifuges
- 3. Thin walls to accelerate termperature transfer
- 4. Clearly visible with either clear or turbid samples
- 5. Available in five colors: clear, pink, green, blue and yellow

# **IV. CAP DESIGN**

- 1. 0.2ml 8-Strip PCR tubes
  - a) Available in flat-top design
  - b) Available in dome-shaped design
  - c) Available in five colors: clear, pink, green, blue and yellow

#### V. PACKAGING

1. 0.2ml 8-Strip PCR tubes: 125 strips/pack, 10 packs/case

## **VI. DESCRIPTION**

Item	Description	Packaging	Sterile
60-0088	0.2ml 8-Strip, Flat-top, Clear	125 strips/pack, 10 packs/case	No
60-0188	0.2ml 8-Strip, Flat-top, Pink	125 strips/pack, 10 packs/case	No
60-0288	0.2ml 8-Strip, Flat-top, Green	125 strips/pack, 10 packs/case	No
60-0388	0.2ml 8-Strip, Flat-top, Blue	125 strips/pack, 10 packs/case	No
60-0488	0.2ml 8-Strip, Flat-top, Yellow	125 strips/pack, 10 packs/case	No

Tel: +86-531-67802668 Fax: +86-531-67803768 www.BiologixGroup.com



Document Number: Document Date: Page BD-T0060 May 2012 2/3

60-0008	0.2ml 8-Strip, Dome cap, Clear	125 strips/pack, 10 packs/case	No
60-0108	0.2ml 8-Strip, Dome cap, Pink	125 strips/pack, 10 packs/case	No
60-0208	0.2ml 8-Strip, Dome cap, Green	125 strips/pack, 10 packs/case	No
60-0308	0.2ml 8-Strip, Dome cap, Blue	125 strips/pack, 10 packs/case	No
60-0408	0.2ml 8-Strip, Dome cap, Yellow	125 strips/pack, 10 packs/case	No

# **VII. STERILITY**

1. N/A

# VIII. CERTIFICATION

Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.


## **IX. RESTRICTIONS**

Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

# X. PRODUCT PARAMETER

# Dimensions=mm

Cat #	Description	Α	В	С	D	E	F	G	Н	I
60-0008	Tube	20.70	3.40	O.D. 8.13	O.D. 7.11	O.D. 6.10	O.D.5.36	9.35	71.12	-
60-0088	Tube	20.70	3.40	O.D. 8.13	O.D. 7.11	O.D. 6.10	O.D.5.36	9.35	71.12	-
60-0008C	Сар	5.51	0.76	O.D. 4.42	O.D. 5.46	80.29	5.51	8.99	7.37	5.52
60-0088C	Сар	3.28	0.76	O.D. 4.42	O.D. 5.46	77.6	3.28	8.99	7.37	5.52




Tel: +86-531-67802668 Fax: +86-531-67803768 www.BiologixGroup.com




Document Number: **BD-T0060** 

Document Date: May 2012

Page 3**/3** 





Fax: +86-531-67803768