

Pipettes

I. APPLICATIONS

1. Widely used in tissue culture as well as in other biological studies.
2. Crucial to the precise handling of liquids.

II. FEATURES

1. Constructed of highly transparent polystyrene
2. Six different volumes, each labeled with a different color for convenience
3. Individually wrapped
4. Non-pyrogenic
5. Non-cytotoxic
6. DNase, RNase, and human DNA free
7. Sterile

III. DESIGN

1. Individually wrapped
2. Color-coded sizes for convenience
3. Ascending and descending graduations

IV. PACKAGING

1. Individually wrapped, 200 Pieces/pack, 800 Pieces/case. (07-5001)
2. Individually wrapped, 150 Pieces/pack, 600 Pieces/case. (07-5002)
3. Individually wrapped, 100 Pieces/pack, 400 Pieces/case. (07-5005, 07-5010)
4. Individually wrapped, 50 Pieces/pack, 200 Pieces/case. (07-5025)
5. Individually wrapped, 40 Pieces/pack, 160 Pieces/case. (07-5050)

V. DESCRIPTION

Cat #	Color	Working Volume (ml)	Packaging
07-5001	Yellow	1.00	200 Pieces/pack, 800 Pieces/case
07-5002	Green	2.00	150 Pieces/pack, 600 Pieces/case
07-5005	Blue	5.00	100 Pieces/pack, 400 Pieces/case
07-5010	Orange	10.00	100 Pieces/pack, 400 Pieces/case
07-5025	Red	25.00	50 Pieces/pack, 200 Pieces/case
07-5050	Purple	50.00	40 Pieces/pack, 160 Pieces/case

VI. STERILITY

1. Sterile

VII. CERTIFICATION


Biologix products are manufactured completely according to the TUV ISO 9001:2008 standard.

VIII. RESTRICTIONS


Refer to the Chemical Resistance Chart in our catalog, and available within our online catalog.

IX. PRODUCT PARAMETER


Dimensions=mm




07-5001




07-5002




07-5005



07-5010



07-5025



07-5050